
Product Data Sheet

Product Name: Fmoc-DL-Ala-OH

Cat. No.: GC65869

Chemical Properties

Cas. No. 35661-38-2

Formula $C_{18}H_{17}NO_4$

M.Wt 311.33

Solubility DMSO : 100 mg/mL (321.20 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-DL-Ala-OH (FMOC-DL-Alanine) is a Fmoc protected alanine derivative^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
