

---

**Product Data Sheet**

---

Product Name: SB1-G-187

Cat. No.: GC65898

**Chemical Properties**

Cas. No.

Formula C<sub>52</sub>H<sub>54</sub>F<sub>3</sub>N<sub>9</sub>O<sub>10</sub>

M.Wt 1022.03

Solubility DMSO : 100 mg/mL (97.84 mM; Need ultrasonic)

Storage Store at -20°C, stored under nitrogen

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**SB1-G-187, a **PROTAC**, is a multi-kinase degrader<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---