

---

**Product Data Sheet**

---

Product Name: TLR7/8 agonist 4

Cat. No.: GC65915

**Chemical Properties**

Cas. No. 2388520-33-8

Formula  $C_{18}H_{24}N_6$  M.Wt 324.42Solubility Storage Store at  $-20^{\circ}C$ 

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

TLR7/8 agonist 4 (compound 41) is a potent **TLR7/8** agonist. TLR7/8 agonist 4 has anti-cancer activity<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---