
Product Data Sheet

Product Name: Thrombin inhibitor 5

Cat. No.: GC65938

Chemical Properties

Cas. No. 328108-09-4

Formula $C_{11}H_9FN_4O_3$ M.Wt 264.21

Solubility DMSO : 31.25 mg/mL (118.28 mM; ultrasonic and warming and heat to 80°C) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

WAY-358871 is a bioactive compound.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA