

---

## Product Data Sheet

---

Product Name: BAY-091  
Cat. No.: GC65952

### Chemical Properties

Cas. No.

Formula  $C_{26}H_{21}FN_4O_2$  M.Wt 440.47

Solubility DMSO : 125 mg/mL (283.79 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

BAY-091 is a novel potent and highly selective inhibitor of the kinase **PIP4K2A**.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA