
Product Data Sheet

Product Name: Sebetralstat

Cat. No.: GC65958

Chemical Properties

Cas. No. 1933514-13-6

Formula C₂₆H₂₆FN₅O₄

M.Wt 491.51

Solubility DMSO : ≥ 125 mg/mL (254.32 mM)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Sebetralstat (KVD900) is a **plasma kallikrein** inhibitor (WO2016083820). Sebetralstat can be used for the research of metabolic diseases^[1].

Sebetralstat is a plasma kallikrein inhibitor^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA