
Product Data Sheet

Product Name: Adenosine receptor antagonist 3

Cat. No.: GC66007

Chemical Properties

Cas. No. 2400864-80-2

Formula $C_{16}H_{14}N_6S$

M.Wt 322.39

Solubility DMSO : ≥ 125 mg/mL (387.73 mM)

Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Adenosine receptor antagonist 3 is a potent antagonist of **adenosine receptor**. Adenosine receptor antagonist 3 has the potential for the research of cancer disease (extracted from patent WO2019233994A1, compound 1)^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA