
Product Data Sheet

Product Name: (S)-1-Benzylpyrrolidin-3-ol

Cat. No.: GC66076

Chemical Properties

Cas. No. 101385-90-4

Formula C₁₁H₁₅NO

M.Wt 177.25

Solubility DMSO : 250 mg/mL (1410.44 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

(S)-1-Benzylpyrrolidin-3-ol (Compound 14) can be used to synthesis impurities of the anti-hypertensive drug Barnidipine hydrochloride [1].

\$

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
