

---

**Product Data Sheet**

---

Product Name: (R)-2-Amino-3-(m-tolyl)propanoic acid

Cat. No.: GC66078

**Chemical Properties**

Cas. No. 114926-39-5

Formula  $C_{10}H_{13}NO_2$  M.Wt 179.22Solubility Storage Store at  $-20^{\circ}C$ 

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

(R)-2-Amino-3-(m-tolyl)propanoic acid is a phenylalanine derivative<sup>[1]</sup>.

Amino acids and amino acid derivatives have been commercially used as ergogenic supplements. They influence the secretion of anabolic hormones, supply of fuel during exercise, mental performance during stress related tasks and prevent exercise induced muscle damage. They are recognized to be beneficial as ergogenic dietary substances<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA