

---

**Product Data Sheet**

---

Product Name: Fmoc-Val-Ala-PAB-PNP

Cat. No.: GC66110

**Chemical Properties**

Cas. No. 1394238-92-6

Formula C<sub>37</sub>H<sub>36</sub>N<sub>4</sub>O<sub>9</sub> M.Wt 680.7

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Val-Ala-PAB-PNP is a cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs)<sup>[1]</sup>.

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---