

---

**Product Data Sheet**

---

Product Name: m-PEG5-CH<sub>2</sub>CH<sub>2</sub>COOH

Cat. No.: GC66152

**Chemical Properties**

Cas. No. 1347750-72-4

Formula C<sub>14</sub>H<sub>28</sub>O<sub>8</sub> M.Wt 324.37

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

m-PEG5-CH<sub>2</sub>CH<sub>2</sub>COOH is a PEG-based based PROTAC linker can be used in the synthesis of PROTACs.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---