
Product Data Sheet

Product Name: Anti- α -Tubulin Antibody, AF555 conjugate

Cat. No.: GC66198

Chemical Properties

Cas. No.

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Anti- α -Tubulin Antibody, AF555 conjugate is a conjugate of mouse anti- α -tubulin monoclonal antibody and the red fluorescent dye Alexa Fluor 555. Anti- α -Tubulin Antibody, AF555 conjugate can be used for the detection of **tubulin** (Ex/Em: 554/567 nm)[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA