
Product Data Sheet

Product Name: Di(N-succinimidyl)adipate

Cat. No.: GC66205

Chemical Properties

Cas. No. 59156-70-6

Formula $C_{14}H_{16}N_2O_8$ M.Wt 340.29

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Di(N-succinimidyl) adipate (DSAP) is an alkyl chain-based PROTAC linker that is applicable to the synthesis of PROTACs.

[1] Maglinao M, et al. J Control Release. 2014 Feb 10;175:36-42.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA