
Product Data Sheet

Product Name: NH-bis(C1-Boc)

Cat. No.: GC66235

Chemical Properties

Cas. No. 85916-13-8

Formula $C_{12}H_{23}NO_4$ M.Wt 245.32

Solubility Storage 4°C, protect from light

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

NH-bis(C1-Boc) is a uncleavable linker used for antibody-drug conjugates (**ADC**).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA