
Product Data Sheet

Product Name: Fmoc-PEG4-Ala-Ala-Asn-PAB

Cat. No.: GC66274

Chemical Properties

Cas. No. 2055048-57-0

Formula $C_{43}H_{55}N_5O_{13}$ M.Wt 849.92

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Fmoc-PEG4-Ala-Ala-Asn-PAB is a cleavable 4 unit PEG ADC linker used in the synthesis of antibody-drug conjugates (ADCs)^[1].

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA