

---

**Product Data Sheet**

---

Product Name: Linagliptin-d4

Cat. No.: GC66355

**Chemical Properties**

Cas. No. 2140263-92-7

Formula C<sub>25</sub>H<sub>24</sub>D<sub>4</sub>N<sub>8</sub>O<sub>2</sub>

M.Wt 476.57

Solubility DMSO : 25 mg/mL (52.46 mM; ultrasonic and warming and heat to 80°C)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Linagliptin-d4 is deuterium labeled Linagliptin. Linagliptin is a highly potent, selective DPP-4 inhibitor with IC<sub>50</sub> of 1 nM.

Stable heavy isotopes of hydrogen, carbon, and other elements have been incorporated into drug molecules, largely as tracers for quantitation during the drug development process. Deuteration has gained attention because of its potential to affect the pharmacokinetic and metabolic profiles of drugs<sup>[1]</sup>.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA