
Product Data Sheet

Product Name: DB-0646
Cat. No.: GC66361

Chemical Properties

Cas. No.

Formula $C_{52}H_{51}F_3N_{14}O_7$

M.Wt 1041.05

Solubility DMSO : 50 mg/mL (48.03 mM; Need ultrasonic)

Storage 4°C, stored under nitrogen

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

DB-0646, a **PROTAC**, is a multi-kinase degrader^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA