
Product Data Sheet

Product Name: Panitumumab

Cat. No.: GC66405

Chemical Properties

Cas. No. 339177-26-3

Formula M.Wt 147kDa

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Panitumumab (anti-EGFR, ABX-EGF) is a recombinant, fully human IgG2 monoclonal antibody that binds to the epidermal growth factor receptor (EGFR).

[1] Wells A Messersmith, et al. Clin Cancer Res . 2007 Aug 15;13(16):4664-6.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA