
Product Data Sheet

Product Name: Etokimab

Cat. No.: GC66407

Chemical Properties

Cas. No. 2022981-44-6

Formula M.Wt

Solubility Storage Store at -80°C. Aliquots should be stored at the same temperature after first use to avoid multiple freeze-thaws

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Etokimab (Antibody ANB 020) is a humanized monoclonal antibody that targets IL-33. Etokimab can be used for the research of atopic dermatitis^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA