
Product Data Sheet

Product Name: CDK9-IN-10

Cat. No.: GC66475

Chemical Properties

Cas. No. 3542-63-0

Formula C₂₂H₁₆O₅

M.Wt 360.36

Solubility DMSO : ≥ 100 mg/mL (277.50 mM)

Storage 4°C, away from moisture and light

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

CDK9-IN-10 is a potent **CDK9** inhibitor. CDK9-IN-10 is the ligand for the PROTAC CDK9 degrader-2 [1].

\$

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA