
Product Data Sheet

Product Name: 4-Amino-L-phenylalanine hydrochloride

Cat. No.: GC66484

Chemical Properties

Cas. No. 62040-55-5

Formula $C_9H_{13}ClN_2O_2$ M.Wt 216.66

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

4-Amino-L-phenylalanine (H-Phe(4-NH₂)-OH) hydrochloride is an endogenous metabolite.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA