
Product Data Sheet

Product Name: Sulfo-SIAB

Cat. No.: GC66617

Chemical Properties

Cas. No. 106145-13-5

Formula $C_{13}H_{11}IN_2O_8S$

M.Wt

482.2

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Sulfo-SIAB is a noncleavable monovalent bilinker.

[1]. Gubala V, et al. A comparison of mono and multivalent linkers and their effect on the colloidal stability of nanoparticle and immunoassays performance. *Talanta*. 2010 Jun 15;81(4-5):1833-9.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA