
Product Data Sheet

Product Name: ((3-Chlorophenyl)sulfonyl)phenylalanine

Cat. No.: GC66669

Chemical Properties

Cas. No. 1449132-27-7

Formula C₁₅H₁₄ClNO₄S

M.Wt

339.79

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

((3-Chlorophenyl)sulfonyl)phenylalanine is a phenylalanine analogue, contains both an amino group and a carboxyl group in its molecule.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA