

---

**Product Data Sheet**

---

Product Name: Mal-PEG2-alcohol

Cat. No.: GC66701

**Chemical Properties**

Cas. No. 34321-81-8

Formula  $C_8H_{11}NO_4$  M.Wt 185.18Solubility Storage Store at  $-20^{\circ}C$ 

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Mal-PEG2-alcohol is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

[1]. An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---