
Product Data Sheet

Product Name: Biotin-PEG3-NHS ester

Cat. No.: GC66703

Chemical Properties

Cas. No. 1253286-56-4

Formula $C_{23}H_{36}N_4O_9S$ M.Wt 544.62

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Biotin-PEG3-NHS ester is a PEG-based **PROTAC linker** can be used in the synthesis of PROTAC.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein.

PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA