
Product Data Sheet

Product Name: Fmoc-alpha-allyl-L-alanine

Cat. No.: GC66714

Chemical Properties

Cas. No. 288617-71-0

Formula $C_{21}H_{21}NO_4$ M.Wt 351.4

Solubility Storage 4°C, stored under nitrogen

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Fmoc-alpha-allyl-L-alanine is an Fmoc-protected alanine derivative.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA