

---

**Product Data Sheet**

---

Product Name: Fmoc-Phe(3-Me)-OH

Cat. No.: GC66736

**Chemical Properties**

Cas. No. 211637-74-0

Formula C<sub>25</sub>H<sub>23</sub>NO<sub>4</sub> M.Wt 401.45

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-Phe(3-Me)-OH is a Fmoc-protected phenylalanine derivative<sup>[1]</sup>.

[1]. Tanada Mikimasa, et al. Cyclic peptide compound having kans inhibiting action. WO2021090855

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---