
Product Data Sheet

Product Name: N3-PEG3-C2-PFP ester

Cat. No.: GC66829

Chemical Properties

Cas. No. 1807530-07-9

Formula $C_{15}H_{16}F_5N_3O_5$

M.Wt

413.3

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

N3-PEG3-C2-PFP ester is a noncleavable 3-unit PEG linker used in the synthesis of antibody-drug conjugates (ADCs).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA