

---

**Product Data Sheet**

---

Product Name: Acid-C3-SSPy

Cat. No.: GC66883

**Chemical Properties**

Cas. No. 250266-79-6

Formula  $C_9H_{11}NO_2S_2$  M.Wt 229.32

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Acid-C3-SSPy is a cleavable DBA-DM4 linker used in the synthesis of DBA-DM4 , a drug-linker conjugate for ADC<sup>[1]</sup>.

[1]. Samantha Congreve, et al. Antibody drug conjugates (ADC) Current status and mapping of ADC:s in clinical programs.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA

---