
Product Data Sheet

Product Name: m-PEG3-NHS ester

Cat. No.: GC67142

Chemical Properties

Cas. No. 876746-59-7

Formula $C_{12}H_{19}NO_7$ M.Wt 289.28

Solubility Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

m-PEG3-NHS ester is a PEG-based **PROTAC linker** can be used in the synthesis of PROTACs.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

[1]. Lepage ML, et al. Design, synthesis and photochemical properties of the first examples of iminosugar clusters based on fluorescent cores. Beilstein J Org Chem. 2015 May 6;11:659-67.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA