

---

## Product Data Sheet

---

Product Name: m-PEG7-CH<sub>2</sub>CH<sub>2</sub>CHO

Cat. No.: GC67160

### Chemical Properties

Cas. No. 1234369-95-9

Formula C<sub>18</sub>H<sub>36</sub>O<sub>9</sub> M.Wt 396.47

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

### Background

m-PEG7-CH<sub>2</sub>CH<sub>2</sub>CHO is a non-cleavable ADC linker used in the synthesis of antibody-drug conjugates (ADCs). m-PEG7-CH<sub>2</sub>CH<sub>2</sub>CHO is also a PEG-based PROTAC linker that can be used in the synthesis of PROTACs.

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

[1]. Erwin R. Boghaert, et al. Anti-egfr antibody drug conjugates. WO2017214282A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA