
Product Data Sheet

Product Name: m-PEG5-Ms

Cat. No.: GC67171

Chemical Properties

Cas. No. 130955-38-3

Formula $C_{12}H_{26}O_8S$ M.Wt 330.4

Solubility Storage Store at $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at $37^{\circ}C$ and shake it in the ultrasonic bath for a while. Stock solution can be stored below $-20^{\circ}C$ for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

m-PEG5-Ms is a PEG-based **PROTAC linker** can be used in the synthesis of PROTACs. m-PEG5-Ms is a cleavable **ADC linker** used in the synthesis of antibody-drug conjugates (ADCs)^[1].

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

ADCs are comprised of an antibody to which is attached an ADC cytotoxin through an ADC linker.

[1]. Sanxing Sun, et al. Triazolotriazine derivatives as a2a receptor antagonists. WO2020002969A1.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA