
Product Data Sheet

Product Name: m-PEG8-Tos

Cat. No.: GC67178

Chemical Properties

Cas. No. 82217-01-4

Formula $C_{24}H_{42}O_{11}S$ M.Wt 538.65

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure **Background**

Tos-PEG8-m is a derivative of silybin ethers, extracted from patent CN105037337A (compound III-d). Tos-PEG8-m is a PEG-based **PROTAC linker** can be used in the synthesis of Silymarin [1].

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein.

PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins.

[1]. Yuxia Zhao, et al. Derivative of silybin ethers, and synthetic method and application thereof. Patent CN105037337A.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA