

---

## Product Data Sheet

---

Product Name: Boc-Phe-(Alloc)Lys-PAB-PNP

Cat. No.: GC67188

### Chemical Properties

Cas. No. 1160844-44-9

Formula  $C_{38}H_{45}N_5O_{11}$  M.Wt 747.79

Solubility Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Boc-Phe-(Alloc)Lys-PAB-PNP is a used as a cleavable linker for antibody-drug conjugates (**ADC**).

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA