
Product Data Sheet

Product Name: NO2-SPP
Cat. No.: GC67189

Chemical Properties

Cas. No. 663598-61-6

Formula $C_{14}H_{15}N_3O_6S_2$ M.Wt 385.42

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

NO2-SPP is a cleavable linker that is used for making antibody-drug conjugate (**ADC**).

[1]. Qiming Chen, et al. Anti-integrin immunoconjugates, methods and uses. WO2006062779A2.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA