

---

## Product Data Sheet

---

Product Name: Aminoxy-PEG1-azide

Cat. No.: GC67240

### Chemical Properties

Cas. No. 2100306-70-3

Formula  $C_4H_{10}N_4O_2$  M.Wt 146.15

Solubility Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution: ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

### Background

Aminoxy-PEG1-azide is a PEG-based **PROTAC linker** can be used in the synthesis of PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

[1]. Zhang F, et al. Discovery of a new class of PROTAC BRD4 degraders based on a dihydroquinazolinone derivative and lenalidomide/pomalidomide. *Bioorg Med Chem.* 2020 Jan 1;28(1):115228.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA