
Product Data Sheet

Product Name: 4-N3Pfp-NHS ester

Cat. No.: GC67270

Chemical Properties

Cas. No. 126695-58-7

Formula $C_{11}H_4F_4N_4O_4$ M.Wt 332.17

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

4-N3Pfp-NHS ester is a noncleavable **ADC linker** used in the synthesis of antibody-drug conjugates (ADCs)^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA