

---

## Product Data Sheet

---

Product Name: Tos-PEG1-CH2-Boc

Cat. No.: GC67351

### Chemical Properties

Cas. No. 1643957-24-7

Formula  $C_{15}H_{22}O_6S$  M.Wt 330.4

Solubility Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility, please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution: ship with blue ice. All other available size: ship with RT, or blue ice upon request.

Structure

### Background

Tos-PEG1-CH2-Boc is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

[1]. Snaebjornsson MT, et al. Non-canonical functions of enzymes facilitate cross-talk between cell metabolic and regulatory pathways. *Exp Mol Med.* 2018 Apr 16;50(4):34.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA