
Product Data Sheet

Product Name: Fmoc-D-3-Methylphenylalanine

Cat. No.: GC67355

Chemical Properties

Cas. No. 352351-64-5

Formula $C_{25}H_{23}NO_4$ M.Wt 401.45

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Fmoc-D-3-Methylphenylalanine is an alanine derivative^[1].

[1]. Wei Yang, et al. Diamine-linked receptor-specific cyclic peptides. Patent WO2021158463A1

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA