

---

**Product Data Sheet**

---

Product Name: Sulfo-SIAB sodium

Cat. No.: GC67483

**Chemical Properties**

Cas. No. 144650-93-1

Formula  $C_{13}H_{10}IN_2NaO_8S$ 

M.Wt

504.19

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Sulfo-SIAB sodium is a noncleavable monovalent bilinker.

[1]. Gubala V, et al. A comparison of mono and multivalent linkers and their effect on the colloidal stability of nanoparticle and immunoassays performance. *Talanta*. 2010 Jun 15;81(4-5):1833-9.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA