
Product Data Sheet

Product Name: Thalidomide-NH-PEG4-Ms

Cat. No.: GC67505

Chemical Properties

Cas. No. 2140807-24-3

Formula $C_{22}H_{29}N_3O_{10}S$ M.Wt 527.54

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Thalidomide-NH-PEG4-Ms is an E3 ligase ligand-linker conjugate that incorporates Thalidomide based cereblon ligand and a linker used for PROTAC BCL-XL degrader XZ739^[1].

XZ739, a CRBN-dependent **PROTAC BCL-XL degrader** with a **DC₅₀** value of 2.5 nM in MOLT-4 cells after 16 h treatment. XZ739 also induces cell death through caspase-mediated **apoptosis**^[1].

[1]. Xuan Zhang, et al. Discovery of PROTAC BCL-X L Degraders as Potent Anticancer Agents With Low On-Target Platelet Toxicity. Eur J Med Chem. 2020 Apr 15;192:112186.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA