
Product Data Sheet

Product Name: Thalidomide-PEG4-NH2 hydrochloride

Cat. No.: GC67507

Chemical Properties

Cas. No. 2387510-82-7

Formula C₂₁H₂₈ClN₃O₈

M.Wt

485.92

Solubility

Storage

Store at -20°C

General tips

For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Thalidomide-PEG4-NH2 hydrochloride is a synthesized **E3 ligase ligand-linker conjugate** that incorporates the Thalidomide based cereblon ligand and a linker used in PROTAC technology^[1].

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA