

---

## Product Data Sheet

---

Product Name: m-PEG6-C6-phosphonic acid ethyl ester

Cat. No.: GC67565

### Chemical Properties

Cas. No. 2028281-88-9

Formula  $C_{23}H_{49}O_9P$  M.Wt 500.6

Solubility Storage Store at  $-20^{\circ}C$

General tips For obtaining a higher solubility , please warm the tube at  $37^{\circ}C$  and shake it in the ultrasonic bath for a while. Stock solution can be stored below  $-20^{\circ}C$  for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

m-PEG6-C6-phosphonic acid ethyl ester is a PEG-based PROTAC linker that can be used in the synthesis of PROTACs<sup>[1]</sup>.

PROTACs contain two different ligands connected by a linker; one is a ligand for an E3 ubiquitin ligase and the other is for the target protein. PROTACs exploit the intracellular ubiquitin-proteasome system to selectively degrade target proteins<sup>[1]</sup>.

[1]. An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA