

---

## Product Data Sheet

---

Product Name: FASN-IN-5

Cat. No.: GC67894

### Chemical Properties

Cas. No. 1309805-49-9

Formula  $C_{29}H_{26}N_4O_2$  M.Wt 462.54

Solubility DMSO : 200 mg/mL (432.40 mM; Need ultrasonic) Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

FASN-IN-5 (example 11), a **FASN** inhibitor, can be used for the research of TH17- or CSF1 -mediated disease or disorder such as cancer, immunological disorders, and obesity<sup>[1]</sup>.

[1]. David S. MILLAN, et al. Methods of using fasn inhibitors. WO2017189613A1.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA