
Product Data Sheet

Product Name: Aluminum Glycinate

Cat. No.: GC68137

Chemical Properties

Cas. No. 13682-92-3

Formula $C_2H_6AlNO_4$ M.Wt 135.05

Solubility Storage 4°C, away from moisture

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

Aluminum Glycinate, an organo-metallic compound, is an antacid. Aluminum Glycinate can be used for the research of indigestion, acid reflux and ulcers^[1].

[1]. MORET BE, et, al. [Experimental and clinical evaluation of aluminum glycinate, a new intragastric buffer solution for achievement of prolonged physiological pH value]. Arztl Wochensh. 1955 Aug 26;10(34-35):788-90.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA
