
Product Data Sheet

Product Name: D-Phenylalanine-d5

Cat. No.: GC68142

Chemical Properties

Cas. No. 362049-55-6

Formula C₉H₆D₅NO₂

M.Wt 170.22

Solubility H₂O : ≥ 12.5 mg/mL (73.43 mM)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure **Background**

D-Phenylalanine-d5 is the deuterium labeled D-Phenylalanine. D-Phenylalanine is the synthetic dextro isomer of phenylalanine. D-Phenylalanine inhibits biofilm development of *Pseudoalteromonas* sp. SC2014[1].

Stable heavy isotopes of hydrogen, carbon, and other elements have been incorporated into drug molecules, largely as tracers for quantitation during the drug development process. Deuteration has gained attention because of its potential to affect the pharmacokinetic and metabolic profiles of drugs^[1].

[1]. Russak EM, et al. Impact of Deuterium Substitution on the Pharmacokinetics of Pharmaceuticals. *Ann Pharmacother*. 2019;53(2):211-216.

[2]. Li E, et al. D-Phenylalanine inhibits biofilm development of a marine microbe, *Pseudoalteromonas* sp. SC2014. *FEMS Microbiol Lett*. 2016 Sep;363(18).

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA