
Product Data Sheet

Product Name: Ticlatone
Cat. No.: GC68144

Chemical Properties

Cas. No. 70-10-0

Formula C₇H₄ClNOS

M.Wt 185.63

Solubility DMSO : 125 mg/mL (673.38 mM; ultrasonic and warming and heat to 80°C)

Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Ticlatone is an antifungal that can be used for the research of mycoses^[1].

[1]. Lundell E, et, al. [Clinical studies on landromil (tiklaton) in treatment of mycoses]. Mykosen. 1971 Nov 1;14(11):531-3.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA