
Product Data Sheet

Product Name: Ponegromab

Cat. No.: GC68339

Chemical Properties

Cas. No. 2368950-15-4

Formula M.Wt 146.08 kDa

Solubility Storage -80°C, Aliquots should be stored at the same temperature after first use to avoid multiple freeze-thaws.

General tips For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT, or blue ice upon request.

Structure

Background

Ponegromab (PF 06946860) is a potent and selective humanized anti-**GDF15** antibody inhibitor with anti-cachexia activity. Ponegromab binds to GDF15 and prevents the binding of GDF15 to GFRAL, thereby blocking GDF15/GFRAL-mediated signaling.

Ponegromab can be used in the research of cancers^[1].

[1]. E.J. Roeland, et al. 1696TiP Phase Ib study to assess the effect of PF-06946860 on appetite following subcutaneous administration in patients with anorexia and advanced cancer. VOLUME 32, SUPPLEMENT 5, S1185, SEPTEMBER 01, 2021. Annals of oncology.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA