

---

## Product Data Sheet

---

Product Name: Ieramilimab

Cat. No.: GC68367

### Chemical Properties

Cas. No.

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

### Background

Ieramilimab (LAG525; IMP701) is a humanized IgG4 monoclonal antibody that binds to **LAG-3**, resulting in inhibition of LAG-3 interaction with MHC-II molecules<sup>[1]</sup>.

[1]. Patrick Schöffski, et al. Phase I/II study of the LAG-3 inhibitor Ieramilimab (LAG525) ± anti-PD-1 spartalizumab (PDR001) in patients with advanced malignancies. J Immunother Cancer. 2022 Feb;10(2):e003776.

**Caution: Product has not been fully validated for medical applications. For research use only.**

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA