
Product Data Sheet

Product Name: Magrolimab

Cat. No.: GC68368

Chemical Properties

Cas. No. 2169232-81-7

Formula M.Wt

Solubility Storage Store at -20°C

General tips For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20°C for several months.

Shipping Condition Evaluation sample solution : ship with blue ice All other available size: ship with RT , or blue ice upon request.

Structure

Background

Magrolimab (Hu5F9-G4) is a first-in-class humanized anti-**CD47** IgG4 antibody with an anti-tumor activity^[1].

[1]. Caner Saygin, et al. Current and emerging strategies for management of myelodysplastic syndromes. Blood Rev. 2021 Jul;48:100791.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: (909) 407-4943 Fax: (626) 353-8530 E-mail: tech@glpbio.com

Address: 10292 Central Ave. #205, Montclair, CA, USA